

“Designed for  
Outstanding  
Results.”

**Powerful · Flexible · Durable**

**The HSS-88 is Our Best Seller:**  
A powerful wet trimmer with the highest grinding speed. Its quiet, vibration-free motor makes it stand out, and its durability ensures a long lifespan. Need more power? Now also available in an extra-strong version!

**OUR PROVEN FORMULA FOR SUCCESS:**

We combine technology, experience, expertise, and team spirit.



**Always here for**

Do you want to request a quote, order products, or have questions about your order, delivery, or warranty?

We look forward to your call and are happy to assist you.



[www.wassermann.hamburg](http://www.wassermann.hamburg)



[fb.com/WassermannDentalmaschinen](https://fb.com/WassermannDentalmaschinen)

**HSS-88**  
**HSS-88 Super\***

Also available in an orthodontic version.



Fig.: HSS-88 with Diamond Grinding Disc

**Info-Hotline**

+49 (0)40 730 926-0



SCAN ME

# HSS-88 HSS-88 Super\*

Also Available in Orthodontic Version:  
For bite-oriented trimming  
of dental models

QUALITY MADE IN GERMANY SINCE 1927

## Why Choose a Wassermann Product?

Wassermann products are made from high-quality materials with a deliberate minimal use of plastics. Our devices are renowned for their **practicality, functionality, and durability**. This commitment results in a high level of trust from users in the Wassermann brand.

## Key Features at a Glance:

- Compatible with all Wassermann grinding discs
- High-torque, maintenance-free motor
- Powerful disc and interior rinsing
- Sturdy, oxidation-resistant work table with side openings
- Standard automatic water supply
- Extensive expansion options
- Proven metal construction with a refined stainless steel front ensures stability and durability
- Large work area of 10.5 x 13.5 cm

## Technical Specifications

Power: 1,000 W / \*1,500 W  
Current: 4.6 A / \*6.5 A

Base Dimensions (W x H x D.): 330 x 335 x 410 mm  
Weight, without Grinding Disc: 20.0 kg / \*22.5 kg

Sound Level: ≤ 70 dB(A)  
Speed: 3,000 RPM

### HSS-88 Article Numbers:

- 168998: HSS-88 with Diamond Grinding Disc  
168996: HSS-88 with Dia-Quick Grinding Disc, partially diamond-coated  
168999: HSS-88 with Corundum Grinding Disc  
168991: HSS-88 with Foil Grinding Disc  
168997: HSS-88 without Grinding Disc

### HSS-88 Super Article Numbers:\*

- 168988: HSS-88 Super with Diamond Grinding Disc  
168989: HSS-88 Super with Corundum Grinding Disc  
168987: HSS-88 Super without Grinding Disc

### Individual Components for Orthodontic Equipment:

- 168140: HSS-88 Grinding Device KMS for Begg Technique  
168120: HSS-88 Alignment Laser  
168088: HSS-88 Super Alignment Laser  
164034: Template for Upper Jaw to Trim Lower Jaw Tubing

## Info & Consultation:

Wassermann Dental-Maschinen GmbH  
Rudorffweg 15 - 17  
21031 Hamburg, Germany

Tel: +49 (0)40 730 926-0  
www.wassermann-dental.com  
sales@wassermann-dental.com